**RESISTA<sup>®</sup>**

**IBK**

**ONESEAL**

2025

# company profile

**RESISTA®**

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**RESISTA**, trademark owned by Ing. Carlo Alberto Issoglio & C. S.r.l. , identifies the full range of products designed by the company to meet the needs of dentists and dental technicians.

**RESISTA is known all over the world since 1946, and synonymous of guarantee, balance and continuity.**

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The company has one of the most ip-to-date manufacturing site for the medical devices, located in Omegna, a beautiful town on Orta Lake, Italy.

With commitment, merit and pride, the technical and scientific staff of the R&D department is daily involved to promote innovation and the improvement of our products and processes.

Over the years, the **RESISTA** brand has also become synonymous of fast and qualitative response to market changes and the needs of modern dentistry.

Quality, Price and Service summarize in three simple words the daily perception of the users of **RESISTA** branded products.

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**RESEARCH & DEVELOPMENT** - The R&D Resista Group's program stems from indications from the clinical world combined with our experience in the field of implantable medical devices.

Due to the use of 3D modeling programs and sophisticated computer systems, we simulate the final geometries and the design, developing the rapid prototyping phases of the projects.

**INNOVATIONS** - Resista Group is structured to meet all the needs of the metal-implant with an advanced conception technology, prototyping and final product engineering to arrive at innovative solutions. Validation is carried out in collaboration with accredited Research Institutes following the International Standards Requirements.

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**PROJECT AND PROCESS QUALITY** - ICIM Spa certified **Ing. C. A. Issoglio & C. S.r.l.** in accordance with **UNI EN ISO 9001** and **UNI EN ISO 13485** in compliance with all current regulations relating to the products and services offered. Every employee covering different roles (technical, engineering, administrative, commercial, mechanic operators, etc.), follows guidelines and objectives for a single purpose: continuous product improvement.

**MAIN OBJECTIVE** - Customer satisfaction is our objective. The winning factor is the company's ability to quickly resolve any type of request. The technical support provided by Product Specialists, availability, efficiency and friendliness are the strengths of our structure.

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**BIOMECHANICAL CONCEPTS** - The **Resista Group** method success is the result of several product design studies and the high technology applied in the production. The precision and reliability of the implant has a strong influence on the long term clinical success.

**The manufacturing process** is carried out by high qualified operators, skilled enough to obtain optimal and reproducible results with the use of electronically controlled CNC machinery.


The machinery tolerances, especially in **prosthetic connections**, are our **gold standard** and they are made between **7 - 10  $\mu\text{m}$** .

**MATERIALS** - The implants are made in **Titanium ASTM Gr.4 Cold Worked** (ISO 5832/2), the prosthetic screws and the prosthetic components are made in **Titanium ASTM Gr.5 ELI** (ISO 5832/2).

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
**SURFACE TREATMENTS** - The non-contaminant micro-subtraction process, **DAE** (Double Acid Etching), changes the implant micro roughness (**Ra, Rq**) and the surface texture maximizing the wettability and biomimetic properties. The new surface treatment and decontamination process are patented by **Nobil Bio Ricerche** improving the bio-chemical properties of the implants.

**QUALITY SYSTEM (QS)** - **UNI EN ISO 9001** and **UNI EN ISO 13485** in compliance of Medical Devices Directives. We are qualified in the design and in the production management of dental implants, dental prosthesis, intraligamental anesthesia syringes and abrasive discs.

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
## ✓ Hydrocarbon Deposition Attenuates Osteoblast Activity on Titanium


R. Hayashi, T. Ueno, S. Migita, Y. Tsutsumi, H. Doi, T. Ogawa, T. Hanawa, and N. Wakabayashi - J Dent Res 93(7):698-703, 2014

**GLASS AND TITANIUM PACKAGING** - All **Resista** implant lines come triple-sterile packaged, featuring an inner container made of **Glass and Titanium Gr 4**.

These two materials are considered the gold standard for **preserving biochemical properties**, as they remain atomically unalterable during Gamma Irradiation process.

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
The conical hexagonal connection ensures an increase in mechanical resistance to extra axial load, greater contact surface, greater connection stability between abutment and implant, thanks to a reduced microcirculation of the biological fluids inside.

The **IBK Resista** connection it's made of deep mechanical interface, composed of on each side  $5^\circ$  conical portion, an anti-rotation hexagon and an axial tightening screw.

This combination guarantees perfect stability of the abutment, regardless of the length or implants diameter, and a cold mechanical welding.

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## 5° CONICAL CONNECTION WITH ANTIROTATIONAL HEXAGON


The 5° on the male/female cone walls, causes a so close reduction of the mechanical interface that confers a direct binding (known as "cold welding") between the abutment and the implant.

This effect is fully activated at the end of the prosthetic tightening of the passing screw, that reaches 30 Ncm of torque and guarantees a bacterial seal through the reduction of the interface gap to zero.

In order to remove the abutment it will be necessary to use an alternative removal screw, that needs to be screwed in the original one place: its unique function is to push the abutment from the inside vertically and to remove it from its bound position.

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
## The Micro-Nano Rough DAE surface treatment speed up bone healing processes

1. Removes the manufacturing organic residuals along with Argon Plasma Cleaning
2. Smooths the edge and eliminates the micro-defects
3. Increases surface and wettability (hydrophilicity surface), improving the first fibrin bridges adhesion
4. Increases protein adhesion
5. Maximizes cellular adhesion through the Micro-Nano roughness suitable for actin filaments (phyllopods) anchoring
6. Changes the Titanium surface chemical characteristics improving the biocompatibility and increasing the cellular proliferation and vitality

## Surface topography and chemistry interact with the cellular differentiation processes


The implant surface uses the microtopography as a communication language with the host tissue cells. The Plasma Argon Cold Cleaning with a controlled-environment packaging at a strict protocol procedures, play a key-role in the control of adhered endotoxins (main immunological response player to implant surfaces).

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### VEGF/VEGF-R/RUNX2 Upregulation in Human Periodontal Ligament Stem Cells Seeded on Dual Acid Etched Titanium Disk

Francesca Diomede, Guya Diletta Marconi, Marcos F. X. B. Cavalcanti, Jacopo Pizzicannella, Sante Donato Pierdomenico, Luigia Fonticoli, Adriano Piattelli and Oriana Trubiani

Materials 2020, 13, 706; doi:10.3390/ma13030706

## SURFACE TREATMENT DISTRIBUTION ON IMPLANT'S BODY

All the **Resista's** implants have a surface treatment that is locally differentiated. The logic of surface treatment distribution is dictated by the endo-osseous implant position which requires different characteristics dependin on the muco-gingival component presence around.

1. **Machined Surface:** Brushed and added with TiO<sub>2</sub> Antibacterial Anatase
2. **Nano-rough Surface:** MAE Treatment, nano-rough with average Ra less than one micron
3. **Micro-rough Surface:** 3/5μ: DAE Treatment, micro-rough with average Ra between 3/5μ
4. **Micro-rough Surface:** 4/7μ: DAE Slow Treatment, micro-rough with average Ra between 4/7μ

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### The Bacterial Anti-Adhesive Activity of Double-Etched Titanium (DAE) as a Dental Implant Surface

Morena Petrini, Alessandra Giuliani, Emanuela Di Campli, Silvia Di Lodovico, Giovanna Iezzi, Adriano Piattelli and Simonetta D'Ercole<sup>1</sup>  
International Journal of Molecular Sciences 2020, 21, 8315; doi:10.3390/ijms21218315

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## THE ONESEAL'S IMPLANT LINE STRONG POINTS

- 1. Macro Geometry:** A vast range of implant shapes for different surgical needs.
- 2. Micro Geometry:** Last generation differentiated surface treated (DAE Micro and Nano rough), Certified Decontamination in Argon Cold Plasma Reactor with white clean room packaging.
- 3. Morse 5° Connection:** Intimate mechanical interface that provides a direct "cold weld" fit between abutment and implant, with 3 micron bacterial seal in the gap.
- 4. Prosthetic Parts:** A wide range of prosthetic variants with 7 micron machining tolerances on the hexagonal connection, with emerging curvilinear profiles, platform switching, Titanium Gr5 prosthetic screws and components for angled screwing.
- 5. Surgical Instruments:** Surgical Drills **3-Tech**, maximum efficiency, minimum trauma, smaller friction coefficient, **PRO MSD** technology for counterclockwise bone densifying application long, lasting resistance, perfect visibility of laser marking and surgical millimeters stop.

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
### Various bio-mechanical factors affecting heat generation during osteotomy preparation: A systematic review

Chirag J Chauhan<sup>1</sup>, Darshana N Shah<sup>1</sup>, Foram B Sutaria<sup>1</sup>

Indian J Dent Res. Jan-Feb 2018;29(1):81-92. doi: 10.4103/ijdr.IJDR\_729\_16.

## LOGIC RELATION BETWEEN DRILL AND IMPLANT CORE

The rotary instruments design and shape is based on the logical coherence between the implants size, in the external dimensions, internal core and their penetration properties.

The **IC line Cylindrical Implants** have a tapered apical portion 3mm long. The **IK line Conical Implants** have a tapered apical portion 5mm long.

The drills have respectively 2/3 apical sections, suitable for the implant tunnel preparation, with respect to the implant geometry, maintaining the best bone implant contact (BIC).

Each step drill is equipped with a sharp cutting edge to facilitate penetration, minimizing friction and overheating.

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
### Bone classification: clinical-histomorphometric comparison

Trisi P, Rao W. - Clin Oral Implants Res. 1999 Feb;10(1):1-7. doi: 10.1034/j.1600-0501.1999.100101.x.

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## THE CONDITIONING OF SURGICAL PROTOCOLS

The differentiated implant geometries in the macro/micro shape provide for different applications, both according to the surgical/prosthetic requirements (immediate loading, early loading or postponed loading) and according to the conditions of the bone/anatomical component of the patient who needs rehabilitation. However, surgical maneuvers and alternative protocols are often required to improve implant performance, according to the morphological characteristics of the residual bone structures.

The success of these maneuvers, that is often the result of a customized tools and protocols combination, represents the different response to the problem in relation to the personal operating skills and the presence of tools suitable for their application.

## MECHANICAL OSSEODENSIFICATION IN ANTI-CLOCKWISE ROTATION

The osseodensification using rotary instruments is a recent surgical technique for the implant site preparation that can be associated with different protocols; these protocols can be used in those particular anatomical conditions such as poor bone quality and insufficient vertical / horizontal dimensions.

This non-subtractive osteotomic bone-condensation approach produces an increase in peri-osteotomic bone density, saving the tissue and increasing the primary implant stability.

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


## RESISTA PRO MSD SURGICAL DRILLS - Modular Surgical Drilling


The Modular Surgical Drills with variable geometry are 3-section helicoidal drills with progressive radius rake, that can be used clockwise and anticlockwise.

The diameter of each drill changes according to a constant ratio of 0.4mm (30/34/38/42/46/50): this allows the operator to choose the drill according to the bone quality (over-preparation or under-preparation). All the drills are coated with a DLC (Diamond Like Carbon) synthetic diamond solid lubricant that maximizes performance in terms of mechanical strength and friction reduction.

The **anticlockwise rotation**, reversing the involved forces, generates 3 different effects on the implant tunnel perforation, which can revolutionize the known milling logic.


- 1) Anterior and lateral thrust of the bone removed by the tip, and of the liquids as well, blood and physiological water.
- 2) Ejection and "hammer effect" of the drill that produce an improvement in vertical control.
- 3) Reduction of cutting efficiency to protect sensitive anatomical parts.

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

- ✓ 1) Bone remodeling and bone loss around subcrestal implants: Two distinct entities.  
Spinato S, Lombardi T, Stacchi C.  
Int J Periodontics Restorative Dent. 2023;43(4):410-411.
- ✓ 2) Influence of apico-coronal positioning of tissue-level implants on marginal bone stability during supracrestal tissue height establishment: A multi-center prospective study.  
Spinato S, Bernardello F, Lombardi T, Soardi CM, Messina M, Zaffe D, Stacchi C.  
Clin Implant Dent Relat Res. 2022;24(5):611-620.

## transgingival oneabutment approach

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




Recommended milling burs rotation speed 500 - 800 rpm

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









STANDARD Healing Abutment - Titanium Gr. 5

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ANATOMICAL Healing Abutment - Titanium Gr. 5

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ANATOMICAL Large Healing Abutment - Titanium Gr. 5

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Laboratory Analog - Titanium Gr. 5 - Anodized

**AN 1000** H 12 mm - Universal Shape

 15Ncm

## Transfer per impronta Pick Up - Titanio Gr. 5

Pick Up Impression Coping - Titanium Gr. 5

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Closed Tray Impression Coping - Titanium Gr. 5

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Temporary Glueing Abutment - Titanium Gr. 5


**MT 1000** Temporary Glueing Antirotation Hex - H1mm

**MTR 1000** Temporary Glueing Rotating - H1mm

**MOFL 00** Temporary Glueing **FLAT** - Rotating - H1mm

**VPF00** Vite di Ricambio - Spare Screw

**DMRP120S** Standard Hexagonal Driver Ø 1.20

## Moncone Titanio Diritto - Titanio Gr. 5

Straight Titanium Abutment - Titanium Gr. 5

**MDN 1000** Straight - Antirotation Hex - Narrow Gingival Profile H1.5mm

**VIPA00** Vite di Ricambio - Spare Screw

**DMRP120S** Standard Hexagonal Driver Ø 1.20

## Moncone Titanio Diritto Fresabile - Titanio Gr. 5

Straight Millable Titanium Abutment - Titanium Gr. 5

**MD 1000** Straight - Antirotation Hex - Ø 5.0 Gingival Profile H2mm

**MD 1040** Straight - Antirotation Hex - Ø 5.0 Gingival Profile H4mm

**MD 1050** Straight - Antirotation Hex - Ø 5.0 Gingival Profile H5mm

**VIPA00** Vite di Ricambio - Spare Screw

**DMRP120S** Standard Hexagonal Driver Ø 1.20

## Moncone Titanio Diritto Spallato - Titanio Gr. 5

Straight Shouldered Titanium Abutment - Titanium Gr. 5

**MA 1000** Straight - Antirotation Hex - Ø 5.0 Gingival Profile H2mm

**MA 1040** Straight - Antirotation Hex - Ø 5.0 Gingival Profile H4mm

**VIPA00** Vite di Ricambio - Spare Screw

**DMRP120S** Standard Hexagonal Driver Ø 1.20

# prosthetics parts

32Ncm

**Moncone Titanio Dritto Spallato** - Titanio Gr. 5

**Shoulder Straight Titanium Abutment** - Titanium Gr. 5


**MD 6020** Large Ø 6.0 Gingival Profile - H 2 mm

**MD 6040** Large Ø 6.0 Gingival Profile - H 4 mm

**VIPA00** Vite di Ricambio - Spare Screw

**DMRP120S** Standard Hexagonal Driver Ø 1.20

32Ncm

**Moncone Titanio Angolato 15°** - Titanio Gr. 5

**15° Angled Titanium Abutment** - Titanium Gr. 5


**MA 1500** Ø 5.0 Gingival Profile - H 2 mm

**MA 1540** Ø 5.0 Gingival Profile - H 4 mm

**VIPA00** Vite di Ricambio - Spare Screw

**DMRP120S** Standard Hexagonal Driver Ø 1.20

32Ncm

**Moncone Titanio Angolato 25°** - Titanio Gr. 5

**25° Angled Titanium Abutment** - Titanium Gr. 5

**MA 2500** Ø 5.0 Gingival Profile - H 2 mm

**MA 2540** Ø 5.0 Gingival Profile - H 4 mm

**VIPA00** Vite di Ricambio - Spare Screw

**DMRP120S** Standard Hexagonal Driver Ø 1.20

32Ncm

**Link in Titanio con Calcinabile** - Titanio Gr. 5

**Titanium Link with Castable Abutment** - Titanium Gr. 5

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### Overdenture Ball Abutments - Titanium Gr. 5

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### Overdenture Equator Abutment Kit - Titanium Gr. 5 - PVD TiN

**BPEQ1** H 1 mm - Narrow Gingival Profile

**BPEQ2** H 2 mm - Narrow Gingival Profile

**BPEQ3** H 3 mm - Narrow Gingival Profile

**BPEQ4** H 4 mm - Narrow Gingival Profile

**BPEQ5** H 5 mm - Narrow Gingival Profile

# OT Equator - accessories and spare parts

## Cappette OT Equator

### OT Equator Caps

<b>TFQ 1</b>	Strong Retention
<b>TFQ 2</b>	Standard Retention
<b>TFQ 3</b>	Soft Retention
<b>TFQ 4</b>	Extra Soft Retention
<b>TFQ 5</b>	Dental Laboratory Cap

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### Rotating Bonding Technique Base Link - Titanium Gr. 5

<b>ALEQ</b>	Analog
<b>TISEQ</b>	Transfer
<b>F20EQ</b>	Equator Screw
<b>SCEQ</b>	Castable Equator
<b>TFSB</b>	Smart Box Transfer Black

## Strumenti Equator

### Equator Tools

<b>IEQ</b>	Inserter - Extractor
<b>DMEQ</b>	Manual Driver
<b>DCEQ</b>	Contra-angle Driver

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**LOC Abutment Kit** - Abutment + Metal Box + Caps


- BPLC1**      Kit LOC Abutment H1 mm - Ø4.0
- BPLC2**      Kit LOC Abutment H2 mm - Ø4.0
- BPLC3**      Kit LOC Abutment H3 mm - Ø4.0
- BPLC4**      Kit LOC Abutment H4 mm - Ø4.0
- BPLC5**      Kit LOC Abutment H5 mm - Ø4.0

## Accessori e Ricambi LOC Abutment

**LOC Abutment Accessories and Spare Parts**

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**DMRP120S**    Standard Hexagonal Driver Ø 1.20

**Analogo CAD da Laboratorio con Vite - Titanio Gr. 5**

**Laboratory CAD Analog with Screw - Titanium Gr. 5**

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**Rotating Bonding Technique Base Link - Titanium Gr. 5**

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**Prosthetic screws for angled hole - Titanio Gr. 5**

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**MDT6DS**      Manual Torque Driver for angled hole

**DCT6DS**      RA Contr-angle Torque Driver for angled hole

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Open Links for Gluing for Angled Screwing - Titanium Gr. 5

**IPB TBAI1** Gingival Profile H 1.7mm - Antirotation Hex

**IPB TBAI2** Gingival Profile H 2.5mm - Antirotation Hex

**IPB TBRAI1** Gingival Profile H 1.7mm - Rotating

**IPB TBRAI2** Gingival Profile H 2.5mm - Rotating

**MDT6DS** Manual Torque Driver for angled hole

**DCT6DS** RA Contr-angle Torque Driver for angled hole

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

## Pilastro Toronto Dritto - Titanio Gr. 5

Toronto Straight Abutment - Titanium Gr. 5

- MU 0001** Toronto Straight MUA Abutment H 1mm
- MU 0002** Toronto Straight MUA Abutment H 2mm
- MU 0003** Toronto Straight MUA Abutment H 3mm
- MU 0004** Toronto Straight MUA Abutment H 4mm
- MU 0005** Toronto Straight MUA Abutment H 5mm
- MU 0006** Toronto Straight MUA Abutment H 6mm

**DCDM** [Toronto Straight Abutment Driver](#)

## Pilastro Toronto Angolato 17° - Titanio Gr. 5

17° Toronto Angled Abutment - Titanium Gr. 5

32Ncm 

- MU 1700** Toronto Angled 17° MUA Abutment H3
- MU 1750** Toronto Angled 17° MUA Abutment H5

**VIPA00** Vite di Ricambio - Spare Screw

**DMRP120S** [Standard Hexagonal Driver Ø 1.20](#)

## Pilastro Toronto Angolato 30° e 45° - Titanio Gr. 5

30° and 45° Toronto Angled Abutment - Titanium Gr. 5

32Ncm 

- MU 3000** Toronto Angled 30° MUA Abutment H3
- MU 3050** Toronto Angled 30° MUA Abutment H5
- MU 4500** Toronto Angled 45° MUA Abutment H4.5

**VIPA00** Vite di Ricambio - Spare Screw

**DMRP120S** [Standard Hexagonal Driver Ø 1.20](#)

# toronto prosthetics

## Pilastro Temporaneo Toronto Rotante - Titanio Gr. 5 Rotating Toronto Temporary Abutment - Titanium Gr. 5

**TT 1000** Rotating Temporary

VPMU 00 Vite di Ricambio - Spare Screw  
**DMRP120S** Standard Hexagonal Driver Ø 1.20

## Pilastro Temporaneo Toronto Antirotazione - Titanio Gr. 5 Antirotation Hex Toronto Temporary Abutment - Titanium Gr. 5

**TTA 1000** Antirotation Hex Temporary

VPMU 00 Vite di Ricambio - Spare Screw  
**DMRP120S** Standard Hexagonal Driver Ø 1.20

## Pilastro Toronto Calcinabile Rotante Rotating Castable Toronto Abutment

**TC 1000** Toronto Castable Abutment - Rotating

VPMU 00 Vite di Ricambio - Spare Screw  
**DMRP120S** Standard Hexagonal Driver Ø 1.20

## Pilastro Toronto Calcinabile Antirotazione Antirotation Hex Castable Toronto Abutment

**TCA 1000** Toronto Castable Abutment - Antirotation Hex


VPMU 00 Vite di Ricambio - Spare Screw  
**DMRP120S** Standard Hexagonal Driver Ø 1.20

## Pilastro Toronto in Cromo Cobalto CoCr Toronto Abutment

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## Transfer Toronto Rotante - Titanio Gr. 5

Rotating Toronto Impression Coping - Titanium Gr. 5

**MU 10TR** Rotating Toronto Impression Coping

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
## Analogo CAD da Laboratorio con Vite - Titanio Gr. 5

Laboratory CAD Analog with Screw - Titanium Gr. 5

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Prosthetic screw for angled hole - Titanio Gr. 5

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
**MDT6DS** Manual Torque Driver for angled hole  
**DCT6DS** RA Contr-angle Torque Driver for angled hole

# oneabutment octa prosthetics

32Ncm

## Oneabutment Octa 3.8 - Titanio Gr. 5

Oneabutment Octa 3.8 - Titanium Gr. 5

- IPBO 381 Oneabutment Transgengivale H1.5 mm
- IPBO 382 Oneabutment Transgengivale H2.5 mm
- IPBO 383 Oneabutment Transgengivale H3.5 mm
- IPBO 384 Oneabutment Transgengivale H4.5 mm
- IPBO 385 Oneabutment Transgengivale H5.5 mm

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## Oneabutment Octa 4.8 - Titanio Gr. 5

Oneabutment Octa 4.8 - Titanium Gr. 5

- IPBO 481 Oneabutment Transgengivale H1.5 mm
- IPBO 482 Oneabutment Transgengivale H2.5 mm
- IPBO 483 Oneabutment Transgengivale H3.5 mm
- IPBO 484 Oneabutment Transgengivale H4.5 mm
- IPBO 485 Oneabutment Transgengivale H5.5 mm

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
15Ncm

## Transfer per impronta Pick Up - Titanio Gr. 5

Pick Up Impression Coping - Titanium Gr. 5

- EPIAO 38 Impression Coping Octa Ø 3.8 mm e Ø 4.8 mm

EPVLO

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**DMRP120S**

Standard Hexagonal Driver Ø 1.20

# oneabutment octa prosthetics

## Analogo Digitale Octa per Laboratorio - Titanio Gr. 5

Laboratory Octa Digital Analog - Titanium Gr. 5

**EPBALO 38DS** Analogo Digitale Octa Ø 3.8 mm

**EPBALO 48DS** Analogo Digitale Octa Ø 4.8 mm

## Pilastro Temporaneo Octa - Titanio Gr. 5

Octa Temporary Abutment - Titanium Gr. 5

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
Short Bonding Technique Octa Base Link - Titanium Gr. 5

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


Angled CONE Abutment - Titanium Gr. 5

**IPBCON 22L3** 22.5° Angled CONE Abutment H3 mm **Long 6 mm**

**IPBCON 30L3** 30° Angled CONE Abutment H3 mm **Long 6 mm**

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
CONE Straight Low Abutment 4mm - Titanium Gr. 5

**IPBCON 0S1** Straight CONE Abutment H1.5 mm **Short 4 mm**

VIPA00

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**DMRP120S**

Standard Hexagonal Driver Ø 1.20

32Ncm

## Moncone CONE Diritto Alto 6mm - Titanio Gr. 5

CONE Straight High Abutment 6mm - Titanium Gr. 5

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**CONE Gluing Cap H4mm - Titanium Gr. 5**

**CTCON RS**      CONE Gluing Cap **Short 4 mm** - Rotating

**CTCON AS**      CONE Gluing Cap **Short 4 mm** - Anti-rotation

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
**Scan Abutment and CONE Cap - PEEK**

**IACON S**      Scan Abutment/Peek Cap CONE **Short 4 mm**

**IACON L**      Scan Abutment/Peek Cap CONE **Long 6 mm**

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**Laboratory CONE Digital Analog - Titanium Gr. 5**

**ALCON S**      Digital Analog CONE **Short 4 mm**

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
- ST M13** Manual Tissue Punche - Ø3.0mm
- ST M14** Manual Tissue Punche - Ø4.1mm
- ST M15** Manual Tissue Punche - Ø5.1mm

- ST M23** Contrangle Tissue Punche - Ø3.0mm
- ST M24** Contrangle Tissue Punche - Ø4.1mm
- ST M25** Contrangle Tissue Punche - Ø5.1mm

# calibrated surgical drills

## FCC 20 - Cylindricals

**FCC 2007** Cylindrical Drill - Ø 2.0 mm - Long 7 mm

**FCC 2008** Cylindrical Drill - Ø 2.0 mm - Long 8,5 mm

**FCC 2010** Cylindrical Drill - Ø 2.0 mm - Long 10 mm

**FCC 2011** Cylindrical Drill - Ø 2.0 mm - Long 11,5 mm

**FCC 2013** Cylindrical Drill - Ø 2.0 mm - Long 13 mm

## FKC 34 - Tapered

**FKC 3407** Tapered Drill - Ø 2.9 mm - Long 7 mm

**FKC 3408** Tapered Drill - Ø 2.9 mm - Long 8,5 mm

**FKC 3410** Tapered Drill - Ø 2.9 mm - Long 10 mm

**FKC 3411** Tapered Drill - Ø 2.9 mm - Long 11,5 mm

**FKC 3413** Tapered Drill - Ø 2.9 mm - Long 13 mm

## FKC 38 - Tapered

**FKC 3807** Tapered Drill - Ø 3.2 mm - Long 7 mm

**FKC 3808** Tapered Drill - Ø 3.2 mm - Long 8,5 mm

**FKC 3810** Tapered Drill - Ø 3.2 mm - Long 10 mm

**FKC 3811** Tapered Drill - Ø 3.2 mm - Long 11,5 mm

**FKC 3813** Tapered Drill - Ø 3.2 mm - Long 13 mm

# calibrated surgical drills

## FKC 42 - Tapered

**FKC 4207** Tapered Drill - Ø 3.5 mm - Long 7 mm

**FKC 4208** Tapered Drill - Ø 3.5 mm - Long 8,5 mm

**FKC 4210** Tapered Drill - Ø 3.5 mm - Long 10 mm

**FKC 4211** Tapered Drill - Ø 3.5 mm - Long 11,5 mm

**FKC 4213** Tapered Drill - Ø 3.5 mm - Long 13 mm

## FKC 46 - Tapered

**FKC 4607** Tapered Drill - Ø 3.9 mm - Long 7 mm

**FKC 4608** Tapered Drill - Ø 3.9 mm - Long 8,5 mm

**FKC 4610** Tapered Drill - Ø 3.9 mm - Long 10 mm

**FKC 4611** Tapered Drill - Ø 3.9 mm - Long 11,5 mm

**FKC 4613** Tapered Drill - Ø 3.9 mm - Long 13 mm

## FKC 50 - Tapered

**FKC 5007** Tapered Drill - Ø 4.3 mm - Long 7 mm

**FKC 5008** Tapered Drill - Ø 4.3 mm - Long 8,5 mm

**FKC 5010** Tapered Drill - Ø 4.3 mm - Long 10 mm

**FKC 5011** Tapered Drill - Ø 4.3 mm - Long 11,5 mm

**FKC 5013** Tapered Drill - Ø 4.3 mm - Long 13 mm

# COMPUTER GUIDED SURGERY

## GUIDED SURGERY TOOLS


# computer guided surgery

## MASCHERINA CHIRURGICA - Mascherina on demand multi-impianti

**SURGICAL GUIDE** - Multi-implants on demand surgical guide

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## MODELLI - Modelli Stereolitografati

**MODELS** - Stereolithographed models

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
**IMMEDIATE TEMPORARY PROSTHESIS** - Drilled PMMA and Milled Titanium Framework

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




## FG M - Mucotomes

- FG M41 S** Mucotome Drill - Ø 4.1 mm - Small
- FG M45 L** Mucotome Drill - Ø 4.5 mm - Large
- FG M55 L** Mucotome Drill - Ø 5.3 mm - Large

# guided surgery drills

## FG C - Cortical Drills

- FG C35 S** Cortical Drill - Ø 3.5 mm - Small
- FG C40 S** Cortical Drill - Ø 4.0 mm - Small
- FG C35 L** Cortical Drill - Ø 3.5 mm - Large
- FG C40 L** Cortical Drill - Ø 4.0 mm - Large
- FG C50 L** Cortical Drill - Ø 5.0 mm - Large

## FG 20 S - FIRST DRILLS Small Driver

- FGR 2008 S** First Drill - Ø 2.0 mm - Long 8,5 mm
- FGR 2010 S** First Drill - Ø 2.0 mm - Long 10 mm
- FGR 2011 S** First Drill - Ø 2.0 mm - Long 11,5 mm
- FGR 2013 S** First Drill - Ø 2.0 mm - Long 13 mm
- FGR 2015 S** First Drill - Ø 2.0 mm - Long 15 mm

## FG 20 L - FIRST DRILLS Large Driver

- FGR 2008 L** First Drill - Ø 2.0 mm - Long 8,5 mm
- FGR 2010 L** First Drill - Ø 2.0 mm - Long 10 mm
- FGR 2011 L** First Drill - Ø 2.0 mm - Long 11,5 mm
- FGR 2013 L** First Drill - Ø 2.0 mm - Long 13 mm
- FGR 2015 L** First Drill - Ø 2.0 mm - Long 15 mm

# guided surgery drills

## FG 34 L - 3.4 IMPLANT DRILLS Small Driver

- FGR 3408 S** 3.4 Implant Drill - Ø 2.9 mm - Long 8,5 mm
- FGR 3410 S** 3.4 Implant Drill - Ø 2.9 mm - Long 10 mm
- FGR 3411 S** 3.4 Implant Drill - Ø 2.9 mm - Long 11,5 mm
- FGR 3413 S** 3.4 Implant Drill - Ø 2.9 mm - Long 13 mm
- FGR 3415 S** 3.4 Implant Drill - Ø 2.9 mm - Long 15 mm

## FG 34 L - 3.4 IMPLANT DRILLS Large Driver

- FGR 3408 L** 3.4 Implant Drill - Ø 2.9 mm - Long 8,5 mm
- FGR 3410 L** 3.4 Implant Drill - Ø 2.9 mm - Long 10 mm
- FGR 3411 L** 3.4 Implant Drill - Ø 2.9 mm - Long 11,5 mm
- FGR 3413 L** 3.4 Implant Drill - Ø 2.9 mm - Long 13 mm
- FGR 3415 L** 3.4 Implant Drill - Ø 2.9 mm - Long 15 mm

## FG 38 S - 3.8 IMPLANT DRILLS Small Driver

- FGR 3808 S** 3.8 Implant Drill - Ø 3.3 mm - Long 8,5 mm
- FGR 3810 S** 3.8 Implant Drill - Ø 3.3 mm - Long 10 mm
- FGR 3811 S** 3.8 Implant Drill - Ø 3.3 mm - Long 11,5 mm
- FGR 3813 S** 3.8 Implant Drill - Ø 3.3 mm - Long 13 mm
- FGR 3815 S** 3.8 Implant Drill - Ø 3.3 mm - Long 15 mm

## FG 38 L - 3.8 IMPLANT DRILLS Large Driver

- FGR 3808 L** 3.8 Implant Drill - Ø 3.3 mm - Long 8,5 mm
- FGR 3810 L** 3.8 Implant Drill - Ø 3.3 mm - Long 10 mm
- FGR 3811 L** 3.8 Implant Drill - Ø 3.3 mm - Long 11,5 mm
- FGR 3813 L** 3.8 Implant Drill - Ø 3.3 mm - Long 13 mm
- FGR 3815 L** 3.8 Implant Drill - Ø 3.3 mm - Long 15 mm

# guided surgery drills

## FG 42 S - 4.2 IMPLANT DRILLS Small Driver

- FGR 4208 S** 4.2 Implant Drill - Ø 3.6 mm - Long 8,5 mm
- FGR 4210 S** 4.2 Implant Drill - Ø 3.6 mm - Long 10 mm
- FGR 4211 S** 4.2 Implant Drill - Ø 3.6 mm - Long 11,5 mm
- FGR 4213 S** 4.2 Implant Drill - Ø 3.6 mm - Long 13 mm
- FGR 4215 S** 4.2 Implant Drill - Ø 3.6 mm - Long 15 mm

## FG 42 L - 4.2 IMPLANT DRILLS Large Driver

- FGR 4208 L** 4.2 Implant Drill - Ø 3.6 mm - Long 8,5 mm
- FGR 4210 L** 4.2 Implant Drill - Ø 3.6 mm - Long 10 mm
- FGR 4211 L** 4.2 Implant Drill - Ø 3.6 mm - Long 11,5 mm
- FGR 4213 L** 4.2 Implant Drill - Ø 3.6 mm - Long 13 mm
- FGR 4215 L** 4.2 Implant Drill - Ø 3.6 mm - Long 15 mm

## FG 46 L - 4.6 IMPLANT DRILLS Large Driver

- FGR 4608 L** 4.6 Implant Drill - Ø 3.9 mm - Long 8,5 mm
- FGR 4610 L** 4.6 Implant Drill - Ø 3.9 mm - Long 10 mm
- FGR 4611 L** 4.6 Implant Drill - Ø 3.9 mm - Long 11,5 mm
- FGR 4613 L** 4.6 Implant Drill - Ø 3.9 mm - Long 13 mm
- FGR 4615 L** 4.6 Implant Drill - Ø 3.9 mm - Long 15 mm

## FG 50 L - 5.0 IMPLANT DRILLS Large Driver

- FGR 5008 L** 5.0 Implant Drill - Ø 4.3 mm - Long 8,5 mm
- FGR 5010 L** 5.0 Implant Drill - Ø 4.3 mm - Long 10 mm
- FGR 5011 L** 5.0 Implant Drill - Ø 4.3 mm - Long 11,5 mm
- FGR 5013 L** 5.0 Implant Drill - Ø 4.3 mm - Long 13 mm

# guided surgery tools

## Guided Mounters

**IP4MG 4S**

Mounter Guidata - Internal Hex ● Small Drive 1.5mm

**IP4MG 4L**

Mounter Guidata - Internal Hex ● Large Drive 5.5mm

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

# instruments & tools

## Cacciaviti Manuali - Manual Screw Drivers

		<b>DMRP120XS</b>	<b>Hexagonal 1.20 - Extra Short</b>
		<b>DMRP120S</b>	<b>Hexagonal 1.20 - Short</b>
		<b>DMRP120L</b>	<b>Hexagonal 1.20 - Long</b>
		<b>ACESTR</b>	<b>Estrattore Protesico Prosthetic Extractor</b>
		<b>DM T6DS</b>	<b>T6 Angled Torx Driver - 15mm</b>
		<b>DMEQ</b>	<b>Equator Squared - Short - 15mm</b>

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		<b>ACDRCA</b>	<b>Hexagonal 1.20 - Short</b>
		<b>ACDRCL</b>	<b>Hexagonal 1.20 - Long</b>
		<b>DC T6DS</b>	<b>T6 Angled Torx Driver</b>
		<b>DCEQ</b>	<b>Equator Squared - Medium</b>

# instruments & tools

## Cacciaviti Pilastrì Transgengivali - Transgingival Abutment Drivers

**DCM** Straight Toronto MUA Contra-angle Driver

**MU10AM** Straight Toronto MUA Manual Driver

**DCOC01** **Octagonal** One Abutment Driver - Short

**DCOC02** **Octagonal** One Abutment Driver - Long

**ACOC01** **Octagonal** One Abutment Driver - Short

**ACOC02** **Octagonal** One Abutment Driver - Long

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**ACMOCA** Contra-angle Implant Driver - Short

**ACMOCL** Contra-angle Implant Driver - Long

Hex 2.3mm

**ACMOMC** Manual Implant Driver - Short

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


### Bone Expanders and Compactors - Titanium Gr. 5

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
**KIT 007** Implant Extractor Kit - 2 Extractors + 1 Box + CRND

## Siringhe per Anestesia - Inox

### Anesthesia Syringes - Stainless Steel

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