MBP

MARGINAL BONE
PRESERVATION

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
1. Galindo-Moreno P, León-Cano A, Ortega-Oller I, Monje A, O Valle F, Catena A. Marginal bone loss as success criterion in implant dentistry: Beyond 2 mm. *Clin Oral Implants Res.* 2015;26(4):e28-e34.
2. Galindo-Moreno P, Catena A, Pérez-Sayáns M, Fernández-Barbero JE, O'Valle F, Padial-Molina M. Early marginal bone loss around dental implants to define success in implant dentistry: A retrospective study. *Clin Implant Dent Relat Res.* 2022;24(5):630-642.
3. Windael S, Collaert B, De Buyser S, De Bruyn H, Vervaeke S. Early peri-implant bone loss as a predictor for peri-implantitis: A 10-year prospective cohort study. *Clin Implant Dent Relat Res.* 2021;23(3):298-308.
4. Donnenfeld OW, Marks R, Glickman I. The apically repositioned flap – A clinical study. *J Periodontol.* 1964;35(5):381-387.
5. Tavtigian R. The height of the facial radicular alveolar crest following apically positioned flap operations. *J Periodontol.* 1970;41(7):412-418.
6. Spray JR, Stage, Morris HF, Ochi S. The influence of bone thickness on facial marginal bone response: stage 1 placement through stage 2 uncovering. *Ann Periodontol.* 2000;5(1):119-128.
7. Cicciù M, Pratella U, Fiorillo L, Bernardello F, Perillo F, Rapani A, Stacchi C, Lombardi T. Influence of buccal and palatal bone thickness on post-surgical marginal bone changes around implants placed in posterior maxilla: A multi-centre prospective study. *BMC Oral Health.* 2023;23(1):309.
8. Eriksson R, Albrektsson T. Temperature thresholds for heat induced bone tissue injury: A vitalmicroscopic study in the rabbit. *J Prosthet Dent.* 1983;50(1):101-107.
9. Norton M. Primary stability versus viable constraint: A need to redefine. *Int J Oral Maxillofac Implants.* 2013;28(1):19-21.
10. Cha JY, Pereira MD, Smith AA, Houschyar KS, Yin X, Mouraret S, Brunski JB, Helms JA. Multiscale analyses of the bone-implant interface. *J Dent Res.* 2015;94(3):482-490.
11. Teixeira AB, Shimano AC, Macedo AP, Valente ML, dos Reis AC. Influence of torsional strength on different types of dental implant platforms. *Implant Dent.* 2015;24(3):281-286.
12. Ericsson I, Persson LG, Berglundh T, Marinello CP, Lindhe J, Klinge B. Different types of inflammatory reactions in peri-implant soft tissue. *J Clin Periodontol.* 1995;22(3):255-261.
13. Canullo L, Penarrocha-Oltra D, Soldini C, Mazzocco F, Penarrocha M, Covani U. Microbiological assessment of the implant-abutment interface in different connections. Cross-sectional study after 5 years of functional loading. *Clin Oral Implants Res.* 2015;26(4):426-434.
14. Atieh MA, Ibrahim HM, Atieh AH. Platform switching for marginal bone preservation around dental implants: A systematic review and meta-analysis. *J Periodontol.* 2010;81(10):1350-1366.
15. Kang MH, Jung UW, Cho KS, Lee JS. Retrospective radiographic observational study of 1692 Straumann tissue-level dental implants over 10 years. II. Marginal bone stability. *Clin Implant Dent Relat Res.* 2018;20(5):875-881.
16. Oh TJ, Yoon J, Misch CE, Wang HL. The causes of early implant bone loss: myth or science? *J Periodontol.* 2002;73(3):322-333.
17. Misch CE, Bidez MW, Strong JT. Scientific rationale for dental implant design. In: Misch CE, editor. *Contemporary Implant Dentistry.* 3rd ed. St. Louis: Mosby; 2008. pp. 200-229.
18. Bozkaya D, Muftu S, Muftu A. Evaluation of load transfer characteristics of five different implants in compact bone at different load levels by finite elements analysis. *J Prosthet Dent.* 2004;92(6):523-530.
19. Tomasi C, Tessarolo F, Caola I, Wennström J, Nollo G, Berglundh T. Morphogenesis of peri-implant mucosa revisited: An experimental study in humans. *Clin Oral Implants Res.* 2014;25(9):997-1003.
20. Judgar R, Giro G, Zenobio E, Coelho PG, Feres M, Rodrigues JA, Mangano C, Iezzi G, Piattelli A, Shibli JA. Biological width around one- and two-piece implants retrieved from human jaws. *Biomed Res Int.* 2014;2014:850120.
21. Berglundh T, Lindhe J. Dimension of the periimplant mucosa. Biological width revisited. *J Clin Periodontol.* 1996;23(10):971-973.
22. Linkevičius T, Puisys A, Steigmann M, Vindasiute E, Linkevičienė L. Influence of vertical soft tissue thickness on crestal bone changes around implants with platform switching: A comparative clinical study. *Clin Implant Dent Relat Res.* 2015;17(6):1228-1236.
23. Verardi S, Orsini M, Lombardi T, Ausenda F, Testori T, Pulici A, Oreglia F, Valente NA, Stacchi C. Comparison between two different techniques for peri-implant soft tissue augmentation: Porcine dermal matrix graft versus tenting screw. *J Periodontol.* 2020;91(8):1011-1017.

24. Spinato S, Lombardi T, Stacchi C. Bone remodeling and bone loss around subcrestal implants: Two distinct entities. *Int J Periodontics Restorative Dent.* 2023;43(4):410-411. Spinato S, Bernardello F, Lombardi T, Soardi CM, Messina M, Zaffe D, Stacchi C. Influence of apicocoronal positioning of tissue-level implants on marginal bone stability during supracrestal tissue height establishment: A multi-center prospective study. *Clin Implant Dent Relat Res.* 2022;24(5):611-620.
25. Galindo-Moreno P, León-Cano A, Ortega-Oller I, Monje A, Suárez F, ÓValle F, Spinato S, Catena A. Prosthetic abutment height is a key factor in peri-implant marginal bone loss. *J Dent Res.* 2014;93(7 Suppl):80S-85S.
26. Spinato S, Stacchi C, Lombardi T, Bernardello F, Messina M, Zaffe D. Biological width establishment around dental implants is influenced by abutment height irrespective of vertical mucosal thickness: A cluster randomized controlled trial. *Clin Oral Implants Res.* 2019;30(7):649-659.
27. Becker K, Mihatovic I, Golubovic V, Schwarz F. Impact of abutment materials and dis-/reconnection on soft and hard tissue changes at implants with platform-switching. *J Clin Periodontol.* 2012;39(8):774-780. *Implants Res.* 2019;30(7):649-659.
28. Tallarico M, Caneva M, Meloni SM, Khanari E, Covani U, Canullo L. Definitive abutments placed at implant insertion and never removed: is it an effective approach? A systematic review and metaanalysis of randomized controlled trials. *J Oral Maxillofac Surg.* 2018;76(2):316-324.
29. Vatenas I, Linkevicius T. One abutment one time vs. repeatable abutment disconnections in implants, restored with cemented / screw retained fixed partial dentures: Marginal bone level changes. A systematic review and meta-analysis. *Stomatologija.* 2021;23(2):35-40.
30. Staubli N, Walter C, Schmidt JC, Weiger R, Zitzmann NU. Excess cement and the risk of peri-implant disease: A systematic review. *Clin Oral Implants Res.* 2017;28(10):1278-1290.
31. Wilson TG Jr. The positive relationship between excess cement and peri-implant disease: A prospective clinical endoscopic study. *J Periodontol.* 2009;80(9):1388-1392.
32. Soulami S, Slot DE, van der Weijden F. Implant-abutment emergence angle and profile in relation to peri-implantitis: A systematic review. *Clin Exp Dent Res.* 2022;8(4):795-806. *J Periodontol.* 2009;80(9):1388-1392.

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