(IT)

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(EN)

SURGICAL INSTRUMENTS AND ACCESSORIES

1.PRODUCT IDENTIFICATION

The instruments produced by **Ing. Carlo Alberto Issoglio and C. srl** are medical devices intended for use in the oral cavity with the following functions:

- preparation with cutting of the implant sites in the maxilla and/or mandible, for **Resista ® implants** (also includes all the tools that allow the removal of the already integrated on-site systems)
- preparation with expansion of the implant sites in the maxilla and/or mandible, for **Resista ® implants**
- insertion, tightening and unscrewing of all implants and connecting screws (surgical screws, transmucosal healing screws, through screws for abutments , transfer screws, etc.).

All components of the surgical instruments are sold in non-sterile packaging and are reusable after cleaning and sterilization. Any other use by the doctor is to be considered improper.

2.INTENDED USE

The instruments in question are intended for mechanical use, that is, they have a shank with attachment (UNI EN ISO 1797-1 standard) for contra-angle and must be used with a suitable micromotor. The devices are intended for specialist medical personnel in the dental field. These instruments must be used in specific environments such as dental offices, hospitals and clinics. The operating protocols have been known and consolidated for some time.

3.DESCRPTION

a. Burs/Extensions/Burlier/Tapers/Stops: These are sharp instruments, to be used for the preparation of surgical sites for dental implants.

Their task is to: prepare the implant site hole, tap the site, prepare the implant collar housing flare, remove the implant from the bone, osteoplasty.

Some types have depth notches to allow the determination of the working depth, they have or do not have the predisposition to attach a depth stop. Refer to the catalogues of the individual implant systems for detailed technical specifications.

b. Drivers/Screwdrivers: They are to be used for the insertion, tightening and unscrewing of all implants and connecting screws of the **Resista ® implant lines**. The usage specifications (turns and tightening torques) are indicated in the various surgical protocols reported in the product catalogue/surgical manuals. Some instruments are equipped with an o-ring that allows retention between different instruments and guarantees their transport and use in the oral cavity without the risk of them falling. Since these o-rings are made of plastic material, it is recommended to periodically check their state of wear and sealing capacity. It is recommended to replace them whenever necessary. Worn o-rings can be removed with a simple probe; new ones must be inserted onto the instrument and gently pushed into their seat. Check that they are fully seated and well positioned.

c. Mucotomes: These are instruments used to create a circular incision in the mucosa, of the diameter of the implant to be used.

All surgical instruments and accessories by Ing. Carlo Alberto Issoglio e C. srl are reusable and non-sterile. For more information and a complete list of available products, please refer to the Resista ® product catalogue

4.RAW MATERIAL USED

Contra-angle instruments produced by **Ing. Carlo Alberto Issoglio e C. srl** do not contain materials of human origin, nor of animal origin. It is recommended to check with patients for any allergies to the substances used. Titanium allergy is a very rare event, but possible. Therefore, it is always necessary to check in advance with patients that they do not have allergies to this material either.

The materials used for the production of surgical instruments are the following:

- AISI 630 / 17-4PH steel
- Sadvik Steel Bioline 1RK91
- Kleinox 4197 / AISI 420 steel
- Titanium ASTM F136 GR.5

GENERAL INFORMATION

1. CONTRAINDICATIONS

The use of drills and surgical accessories for the application of an implant-prosthetic project is contraindicated in patients with systemic diseases with compromised immune system, such as AIDS, patients undergoing medical treatment that could compromise the healing of the implant site, patients who do not follow the correct oral hygiene procedures or patients unable to observe the oral hygiene procedures required after the placement of the implants. In addition, excessive use of tobacco increases the risk of complications.

The drills must be used with caution, especially in cases of low bone density. The site must be adequately prepared after careful and scrupulous study of the implant project by the doctor. In any case, the choice of the type of drill to be used is at the discretion and responsibility of the doctor. The surgical instruments must not be used in a manner other than that for which they were designed and described in the various product catalogues, as well as in this document. Any improper use could compromise the correct functioning of the instrument and could cause damage to surgical equipment, as well as compromise the success of the implant project. **Ing. Carlo Alberto Issoglio e C. srl** declines all responsibility.

2. WARNINGS

Resista ® surgical instruments for and with **Resista ®** dental implants and only after having carefully read the information provided to the user.

- a. Surgical instruments MUST be used exclusively by qualified medical personnel familiar with implantology theory and practice.
- b. It is the doctor's responsibility during the treatment planning phase to choose the implant and the type of prosthetic rehabilitation; consequently, it is the doctor's responsibility to choose the type and size of the implant components and related surgical instruments for the preparation of the site and their use.
- c. Use **ONLY Resista ®** surgical instruments for and with **Resista ®** dental implants
- d. Before use, ALWAYS check the integrity of the instrument, the packaging and the blister pack if new. NEVER use the instrument if there are signs of damage. Replace it immediately with a new one.
- e. The washing and sterilization phase must be carried out with a validated method under the care and responsibility of the using Doctor.
- f. Handle all instruments with particular care, so as not to cause accidental damage that could cause incorrect functioning and, during the implant-prosthetic procedure, handle all devices under aseptic conditions and with sterile gloves.
- g. The doctor is responsible for: diagnosis; treatment planning; surgical and prosthetic technique adopted; washing, sterilization and handling of surgical instruments. Errors and inaccuracies in the preparation of the site and positioning of the implant can lead to the failure of the implant-prosthetic treatment.
- h. In case of disposal, the instrument MUST be treated as hospital waste. It is therefore recommended to dispose of the device only after a regular cleaning process, in a controlled manner and in accordance with the applicable regulations in force, whenever the need arises.

i. **Ing. Carlo Alberto Issoglio e C. srl** It also declines any responsibility, expressed or implicit, relating to direct, indirect or any other type of damages that may be connected to or may derive from any errors, of evaluation or professional practice, made in the use of the various **Resista ®** products. The use of dental implants, surgical treatment and prosthetic rehabilitation must in fact take place under the supervision of a doctor, who therefore assumes full responsibility.

j. Surgical instruments produced by **Ing. Carlo Alberto Issoglio e C. srl** are sold in NON-STERILE packaging. Before using them, they must be cleaned, disinfected and sterilized. Failure to comply with this warning may cause infections to the patient. When handling the devices, both during use and during cleaning and sterilization operations, it is recommended to always use surgical gloves for individual protection from bacterial contamination. Failure to comply with these rules may result in cross-infection.

K. Additional Recommendations: The drills, if used properly, produce a hole consistent with the depth of insertion of the implant. The protocols to be used for preparing the site with the various drills are shown in the implant systems catalogue. It is recommended to use (for all types of drills) the rotation speeds always indicated in the product catalogue to avoid the development of bone necrosis. Pressure should never be applied that is so high as to forcefully stop the rotation of the instrument; this could lead to excessive heat build-up in the tissues involved in the cut with consequent possible bone necrosis and ruin both the instrument and the device used (micromotor). In some cases, the instrument itself could also break. It is recommended to work intermittently to avoid overheating and wear of the working part and undue heat build-up in the tissues involved in the cut. The use of an appropriate cooling liquid is recommended. In the absence of adequate irrigation, bone necrosis may occur.

The consumption of the burs depends largely on the type and density of the bone being milled: harder bone leads to greater wear of the instruments. For greater safety and prudence with respect to the wear resistance of the device, it is recommended that the burs be used for no more than 10 work cycles or sooner if the instruments lose their cutting capacity. The 10 recommended cycles represent an average figure. It is recommended to check the maintenance status of the residual cutting capacity after each intervention.

Ing. Carlo Alberto Issoglio e C. srl assumes no responsibility in case of excessive use. The cutters must never be sharpened before use. Never use damaged, bent or worn tools. There are two types of cutter coatings, based on the type of cutter itself: TiN coating (Titanium Nitride) yellow colour and DLC (Diamond Like Carbon) coating black color. Different versions of cutters are also available; please refer to the product catalogue for a detailed description.

A drill extension is available for use when the total length of the instruments is too short due to the presence of adjacent teeth that do not allow the handpiece head to pass through. If you use it, make sure that the drill shank is properly inserted. Incomplete or inappropriate insertion may cause eccentricity of rotation of the drill and could also cause damage to the handpiece or the drill, as well as creating an altered implant site that is not suitable for the type of drill used. There are devices called DRILL STOPS and SPACER BUSHINGS that allow you to limit the working length of a drill to a predetermined depth. The measurements are shown in the catalogues. The drill stops and spacer bushings are also identified by a colour code that determines the range of drill diameters that can be used with that specific stop. For the uses and couplings of the stops or spacer bushings with the relative drills, refer to the product catalogue

3. WASHING / STERILIZATION / STORAGE

Resista ® surgical components for dental implants are sold in NON-STERILE condition. Before use they must be cleaned, disinfected and sterilized. These processes must be carried out prior to intraoral use. Failure to comply with these instructions may lead to the onset of problems for which **Ing. Carlo Alberto Issoglio e C. srl** declines all responsibility.

Washing: immerse the devices in an ultrasonic bath containing a detergent and disinfectant solution, used at the concentration and for the contact times indicated by the manufacturer on the label. After washing, rinse the products with sterile deionized water and dry with a sterile cloth, and if necessary, compressed air.

Sterilization: transfer the medical devices, contained in suitable bags, to a steam autoclave and treat with a standard parameter cycle (121°C for at least 20 minutes for the first use, 134°C for at least 20 minutes for subsequent uses).

Storage: After sterilization, the product must remain in the sterilization pouches. The pouches must be opened only immediately before reuse. The sterilization pouches are normally able to maintain sterility inside them, unless the wrapping is damaged. Do not use the surgical components if the pouches have defects or damage; in this case, proceed with a new sterilization cycle with a new pouch. The storage period of the sterilized products inside the pouches must not exceed that recommended by the pouch manufacturer.

4. REGULATORY REFERENCES

The design and production of surgical components is carried out in compliance with the most up-to-date directives and harmonized standards in force regarding the materials used, the production processes, the information provided and the packaging.

5. DISCLAIMER AND OBLIGATIONS

Ing. Carlo Alberto Issoglio e C. srl declines any responsibility or obligation, direct or indirect:

- for breakage caused by improper use.
- for any damage that may arise from such breakage.
- for any damage caused to the patient and/or following the intervention, as **Ing. Carlo Alberto Issoglio e C. srl** cannot in any way have control over the methods of carrying out the intervention.

You must report any serious incident that occurs in relation to the medical device supplied by us to the Manufacturer, the Notified Body and the Competent Authority of the Member State in which you are located.

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Ing. C. A. Issoglio & C. S.r.l.