Ing. C. A. ISSOGLIO & C. S.r.l.
Via F.lli D. D. 68
28887 OMEGNA(VB) - Italy
+390323282904 Fax +390323282914
info@resista.it - www.resista.it

CE
0425

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(EN) INSTRUCTIONS FOR THE USE OF THE W2 IMPLANT SYSTEM

WARNING: The W2 implant systems can be exclusively used by Dentists, Odontologists and Maxillofacial Surgeons. Esthetic Implant is integral part of a complete system; therefore the use of original instruments and components is recommended.

MATERIALS: The W2 biphasic and UNI-Q-Mua Implant are made of CP Titanium grade 4, (ASTM-F-67). The W2 monophasic implants are made of Titanium grade 5 (ASTM-F-136).

FIELDS OF USE: Implant systems are suitable for any case of prosthesis-functional and aesthetic rehabilitation of both jawbone and mandibular.

The placing of W2 is recommended when the periodontal, functional and physical conditions of patients are favourable.

PACKAGING OF THE W2 IMPLANTS: The system equipped with "Pre-assembled Mounting Device" (Mount) allows professionals to transport and place the desired implant under absolute sterility conditions, thus avoiding direct handling. The implant equipped with Mount is taken out of the sterile package by means of a screw (manual, jack screw or contra-angle), to be directly placed in the surgical site. The Mount is screwed onto the implant by means of a through screw, to be removed when insertion has been completed. The Mount can also be used as impression-transfer

after the second surgical time. The package includes the covering screw. All packages report the main characteristics on the outside: description, lot number, expiry date, reordering code and colour code for simple identification. **ATTENTION:** In order to get the goods replacement, in case the quality does not meet the manufacturer's requirements or in case of a visibly damaged package, it is necessary to inform in advance the distributor by phone within 8 days from the goods receipt; return the product in its original package by filling in regular shipping documents.

STERILITY - DECONTAMINATION: W2-Implants are subjected to a decontamination, which guarantees surface purity free of potential contaminants (manufacturing particles). The packaging is carried out in a Clean Room, and sterilization is carried out by Beta or Gamma-ray irradiation (25 KGV).

RE-STERILIZATION: The autoclave re-sterilization of implants is not allowed. In this case the manufacturer declines all responsibility.

STORAGE: Keep the W2 systems in their original package until they are placed, in a dry place at room temperature. The implant shall not be used beyond the expiry date reported on the package.

PRECAUTIONS AND SIDE EFFECTS: After the operation avoid any stressful physical activity. The side effects of pure titanium are unknown so far. In spite of the high biocompatibility of pure titanium implants, risks of failure are not to be excluded (fibrointegration), which are typically due to technical procedures or the physiological character of biologic processes, not to the product itself. The possible reasons for such failures may be: incorrect surgery techniques and failure to comply with the absolute contra-indications to implantology.

- CONTRA-INDICATIONS:**
 - Lack of proper bone support (min. 1.5 mm perimetric)
 - Patient in poor health conditions
 - Lack of oral hygiene and cooperation by the patient

COMPLICATIONS: To avoid post-surgery complications it is very important to conduct an accurate investigation-study before starting the treatment.

- Damages to the vascular-nervous bundle in the lower arcade, and subsequent paresthesia that leads to the lip insensitivity, resulting in the need to remove the implant.
- Involvement of the nasal sinus and the nasal-palatine ducts, resulting in sinusitis, peri-implant inflammation due to rhinitides, anesthesia and stiffening of the upper lip.

CLINICAL PREPARATION: Case study, detailed case-history, clinical and radiographic complete analyses.

OPERATING TECHNIQUE: The operating technique is described in the surgical protocol. A high degree of accuracy is necessary in the site preparation through the various steps. The thermal shock due to the friction of the cutters shall be minimized by using successive milling cutters, appropriate number of revolutions, irrigation by sterile pre-cooled solutions and proper recovery times. The taking of the implant from the package has to be carried out in compliance with the asepsis regulations. The implant shall not be used both of W2 brand and compatible tools.

DOCUMENTATION: It is necessary to apply, on the patient's clinical record, the label that comes with the implant, reporting the main characteristics of the inserted screw. Moreover it is recommended to prepare sufficient radiological, clinical and photographic documentation, as well as to fill in the informed consent signed by the patient.

Any serious incident occurred in connection with this medical device must be communicated to the Manufacturer, to the Notified Body and to the Competent Authority of the Member State where you are located.

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