Ing. C. A. Issoglio & C. S.r.l.

Via F.lli Di Dio, 68
28887 Omegna (VB)
Tel. +39.0323.828.004 Fax
+39.0323.828914
info@resista.it

IT) COMPONENTISTICA PROTESICA

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(EN) PROSTHETIC PARTS

WARNING - The prosthetic parts of W2 implant line can only be used by Dentists, Maxillofacial Surgeons and Dental Technicians. Those are part of a complete system, so we recommend the use of original prosthetic parts.

WARNING - The device is disposable. It has not to be reused. The manufacturer assumes no responsibility for any alteration to the biocompatibility's properties of the device after re-sterilization.

Ing. C. A. Issoglio & C. S.r.l.

Via F.lli Di Dio, 68
28887 Omegna (VB)
Tel. +39.0323.828.004 Fax
+39.0323.828914
info@resista.it

MATERIALS: The prosthetic parts are made of Titanium Grade 5 - Golden League - Castable Acrylic - PEEK - Zirconia

FIELDS OF USE: The prosthetic components is indicated in all cases of prosthetic-functional and aesthetic rehabilitation of both the maxillar and the mandibular. The prosthesis is made by single crowns, bridges, full or partial dentures attached to plants by appropriate prosthesis.

TECHNICAL OPERATIONS: The operational techniques for use of healing abutments are explained in the Surgical Protocol (for information on how to receive phone or send an email at the numbers listed above). The indications on the remaining prosthetic part are shown on the W2 catalogue.

PRECAUTIONS: Before the replacement follow the post-operative procedures and the healing time recommended in the Surgical Protocol.

NAMES OF PROSTHETIC PARTS: identification data of prosthetic parts used must appear on the patient's medical record. We also suggest to make, radiological and clinical documentation, as also the complete statistics.

INSTRUCTIONS FOR FIRST USE: The medical devices should undergo the following procedure prior to first use.

- 1) Washing and disinfection: treat the medical devices in ultrasonic tank with a suitable detergent – disinfectant solution prepared and used as indicated by producer.
- 2) Sterilization: transfer the medical devices, contained in suitable pouches, to a steam autoclave and process with a cycle (pre-vacuum) at 121°C for 18 minutes.

WASHING: immerse the devices in an ultrasonic tank containing a detergent and disinfectant solution (for example enzymatic based with quaternary ammonium salts or chlorhexidine) or a specific alkaline detergent, used at the concentration and for the contact times indicated by the manufacturer on the label.

After washing, rinse the products with sterile deionized water and dry with a clean cloth that does not release particles.

STERILIZATION: transfer the medical devices, contained in suitable bags, into a steam autoclave and treat with a cycle with standard parameters (minimum 121°C for at least 15 minutes).

PACKAGING: All items are placed in a heat-sealed polyethylene packaging. The packaging is labeled with the reference product code and lot number. The product is supplied non-sterile.

STORAGE: Store the device in its original packaging until ready to use, in a dry place at room temperature.

You must report any serious incident that occurs in relation to the medical device supplied by us to the Manufacturer, the Notified Body and the Competent Authority of the Member State in which you are located.

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