



INTERNAL TRI-LOBE IR

NARROW - **REGULAR** - **WIDE**

TAPERED
SHAPE



IR



IMPIANTO
SENZA MOUNTER
MOUNTERLESS

INTERNAL TRI-LOBE IR



IR 35



IR 43



IR 50

CODE	IMPLANT	LENGTH	PLATFORM	APEX
IR 3508	Ø 3.5	8.5 mm	Ø 3.5	Ø 2.2
IR 3510	Ø 3.5	10 mm	Ø 3.5	Ø 2.2
IR 3511	Ø 3.5	11,5 mm	Ø 3.5	Ø 2.2
IR 3513	Ø 3.5	13 mm	Ø 3.5	Ø 2.2
IR 4308	Ø 4.3	8,5 mm	Ø 4.3	Ø 2.8
IR 4310	Ø 4.3	10 mm	Ø 4.3	Ø 2.8
IR 4311	Ø 4.3	11,5 mm	Ø 4.3	Ø 2.8
IR 4313	Ø 4.3	13 mm	Ø 4.3	Ø 2.8
IR 4316	Ø 4.3	16 mm	Ø 4.3	Ø 2.8
IR 5008	Ø 5.0	8,5 mm	Ø 5.0	Ø 3.8
IR 5010	Ø 5.0	10 mm	Ø 5.0	Ø 3.8
IR 5011	Ø 5.0	11,5 mm	Ø 5.0	Ø 3.8
IR 5013	Ø 5.0	13 mm	Ø 5.0	Ø 3.8
IR 5013	Ø 5.0	13 mm	Ø 5.0	Ø 3.8