

UNI-Q-MUA prosthetics

Transfer e Anologo UNI-Q-MUA - Titanio Gr. 5

UNI-Q-MUA Impression Coping and Analog - Titanium Gr. 5

EP MAL-Q Analog H 10 mm - Universal Shape

EP MI-Q Rotating Transfer

EP MIA-Q Rotating Transfer

EP MVL-Q Vite di Ricambio Q - Q Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



Pilastro Calcinabile UNI-Q-MUA - Acrilico

UNI-Q-MUA Castable Abutment - Acrylic

EP MC-Q Rotating - Standard Gingival Profile

EP MCA-Q Antirotation Hex - Standard Gingival Profile

EP MV-Q Vite di Ricambio Q - Q Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



Pilastro Temporaneo UNI-Q-MUA - Titanio Gr. 5

UNI-Q-MUA Temporary Abutment - Titanium Gr. 5

EP MT-Q Rotating - Bonding Technique - Incollaggio Passivo

EP MTS-Q Rotating - Welding Technique - Saldatura Intraorale

EP MV-Q Vite di Ricambio Q - Q Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



Pilastro Base CoCr Antirotazione - CoCr

Antirotational CoCr Base Abutment - CoCr

EP MCCA-Q Antirotation Hex - Standard Gingival Profile

Intervallo di Fusione 1350°-1430°
Espansione Termica 13,3 - 13,6 µm/(m°C) (20°-200°)

EP MV-Q Vite di Ricambio Q - Q Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



Componenti Digitali UNI-Q-MUA - Titanio Gr. 5

UNI-Q-MUA Digital Components - Titanium Gr. 5

EP MSAQ-DS Antirotational Scan Abutment

EP MTAQ-DS Antirotational Bonding Technique Base Link

EP MTRQ-DS Rotating Bonding Technique Base Link

EP MV-Q Vite di Ricambio Q - Q Spare Screw



DMX 1/2 - Toronto Abutment Torx Driver



UNI-Q-MUA surgical tools



KIT Q
Kit Base UNI-Q-MUA
BOX 02 + 5 Frese + 12 Accessori



Pini di Parallelismo Chirurgici - Titanio Gr. 5

Paralleling Pins - Titanium Gr. 5

| | |
|--------------|-------------------------------------|
| PIN | Double diameter - 1.8mm / 3mm - 0° |
| PIN17 | Double diameter - 1.8mm / 3mm - 17° |
| PIN32 | Double diameter - 1.8mm / 3mm - 32° |



Pini di Parallelismo Protesici - Titanio Gr. 5

Prosthetic Parallelism Pins - Titanium Gr. 5

| | |
|--------------|----------------------------------|
| PIN-Q | Peek Pin for Parallelism Control |
|--------------|----------------------------------|

20Ncm 



Pilastro di Guarigione UNI-Q-MUA - PE

UNI-Q-MUA Healing Abutment - PE

| | |
|-----------------|------------------------------------|
| EP MHA-Q | H 4 mm - Standard Gingival Profile |
|-----------------|------------------------------------|

| | |
|----------------|------------------------------|
| DMX 1/2 | Toronto Abutment Torx Driver |
|----------------|------------------------------|