


# SPARE SCREWS & MOUNTERS

## REPLACEMENT PRODUCTS

### SCREWS & MOUNTERS

























# spare screws

## ESAGONO INTERNO IK / IA / IC - INTERNAL HEXAGON IK / IA / IC

IP VIT20	H 20mm Transfer Screw - Ø 1.8mm	25Ncm		
IP VIT30	H 30mm Transfer Screw - Ø 1.8mm	25Ncm		
IP VIT	H 8mm Prosthetic Screw - Ø 1.8mm	32Ncm		
IP VITG	H 8mm Prosthetic Screw - Ø 1.8mm - Gold PVD TiN	32Ncm		
IP VITT6-DS	Prosthetic Screw T6 per Avvitamento Angolato - Ø 1.8mm	32Ncm		
IP VITMA	Angled 17°/30° MUA Abutment Screw - Ø 1.8mm	32Ncm		

## ESAGONO ESTERNO ES / EC / EK - EXTERNAL HEXAGON ES / EC / EK

EP VIT20	H 20mm Transfer Screw - Ø 2.0mm	25Ncm		
EP VIT30	H 30mm Transfer Screw - Ø 2.0mm	25Ncm		
EP VIT	H 8mm Prosthetic Screw - Ø 2.0mm	32Ncm		
EP VITG	H 8mm Prosthetic Screw - Ø 2.0mm- Gold PVD TiN	32Ncm		
EP VITT6-DS	Prosthetic Screw T6 per Avvitamento Angolato - Ø 2.0mm	32Ncm		
EP VITMA	Angled 17°/30° MUA Abutment Screw - Ø 2.0mm	32Ncm		
EP MV	Toronto Prosthetic Spare Screw - Ø 1.4mm	15Ncm		
EP MVL	Toronto Impression Spare Screw - Ø 1.4mm	15Ncm		
EP VITS	H 8mm CAD CAM Screw - Ø 2mm Torque Gold	32Ncm		
EP VITB	H 8mm - Ø 2mm External Hexagon - Torque Black DLC	32Ncm		
SP VIT	H 8mm Prosthetic Screw - Ø 2.0mm	32Ncm		

# spare screws

## PLATINUM CONICAL CONNECTION IK-P / IA-P / ES-P

### **NARROW** Connessione Conica - Titanio Gr. 5

Conical Connection **NARROW** - Titanium Gr. 5

**IP VIT3L-P** H 20mm Transfer Screw - Ø 1.6mm

25Ncm  

**IP VIT3-P** H 8mm Prosthetic Screw - Ø 1.6mm

25Ncm  

**DMX 1/2** Toronto Abutment Torx Driver

### **REGULAR** Connessione Conica - Titanio Gr. 5

Conical Connection **NARROW** - Titanium Gr. 5

**IP VIT4L-P** H 20mm Transfer Screw - Ø 2.0mm

32Ncm  

**IP VIT4-P** H 8mm Prosthetic Screw - Ø 2.0mm

32Ncm  

**DMX 1/2** Toronto Abutment Torx Driver

## RICAMBI DI MOUNTERS IMPIANTI - IMPLANT MOUNTERS SPARE PARTS

### **MTA Mounter Transfer Abutment** - Titanio Gr. 5

MTA Mounter Transfer Abutment - Titanium Gr. 5

**EPM 40** H 12mm - Ø 4.7mm External Hexagon

**EPM 50** H 12mm - Ø 5.7mm External Hexagon

**IPM 35** H 12mm - Ø 4.5mm Internal Hexagon

**IPM 40** H 12mm - Ø 4.7mm Internal Hexagon

**IPM 50** H 12mm - Ø 5.5mm Internal Hexagon

32Ncm



32Ncm



**DM 1/2** Standard Hexagonal Driver